

OCPA

INFORMATION REPORT INFORMATION REPORT

CENTRAL INTELLIGENCE AGENCY

This material contains information affecting the National Defense of the United States within the meaning of the Espionage Laws, Title 18, U.S.C. Secs. 793 and 794, the transmission or revelation of which in any manner to an unauthorized person is prohibited by law.

S-E-C-R-E-T
NOFORN

25X1

COUNTRY USSR

REPORT

SUBJECT Model T G-3 Crucible-Type Furnace

DATE DISTR. 28 October 1958

(complete description
of equipment and its
operation)

NO. PAGES 1

REFERENCES

25X1

DATE OF
INFO.PLACE &
DATE ACQ.

SING COPY

SOURCE EVALUATIONS ARE DEFINITIVE. APPRAISAL OF CONTENT IS TENTATIVE.

25X1

Model T G-3 crucible-type electric furnace.

the following information about the furnace: certification; application; brief description and photograph; schematic diagram of furnace connection to supply; specification; preparation of furnace for operation; operation; storage; delivery volume; guarantee; and a list and index of spare parts. Classification of the document when separated from this report is OFFICIAL USE ONLY.

28 JAN 1959

26 JAN 1959

25X1

13

S-E-C-R-E-T
NOFORN

25X1

STATE	X	ARMY	X	NAVY	X	AIR	X	FBI		AEC						
(Note: Washington distribution indicated by "X"; Field distribution by "#".)																

INFORMATION REPORT INFORMATION REPORT

FOR OFFICIAL USE ONLY

МОСГОРСОВНАРХОЗ

УПРАВЛЕНИЕ РАДИОТЕХНИЧЕСКОЙ ПРОМЫШЛЕННОСТИ И ПРИБОРОСТРОЕНИЯ

**ГОСУДАРСТВЕННЫЙ ЗАВОД
„ПЛАТИНОПРИБОР“**


МОСКВА

FOR OFFICIAL USE ONLY

FOR OFFICIAL USE ONLY

PLANNING STATE WORKS

25X1



U.S. AIR FORCE

CONTROL APPARATUS MAKING INDUSTRY

MINISTRY OF INSTRUMENT AND AUTOMATIC

U.S.S.R.

22051/5

FOR OFFICIAL USE ONLY

MODEL T P - 3 CRUCIBLE TYPE ELECTRIC FURNACE

C E R T I F I C A T E

OF WORKS CONTROL DEPARTMENT

Serial No. 4

Supply 220 Volts

Consumption 1850 Watts

Maximum working temperature 1200° C.

Dimensions of working room:

Diameter 72 mm.

Height 175 mm.

Time required to raise temperature
up to 1200° C 120 minutes

Uniformity of temperature within the
furnace working room $\pm 5^{\circ}\text{C}$

Stretch of uniform temperature zone 80 mm.

25X1


- 2 -

MODEL IT-3 CRUCIBLE TYPE ELECTRIC FURNACE
OPERATION AND MAINTENANCE INSTRUCTIONS

1. APPLICATION

The furnace is designed for calcination of sediments and for thermal analysis of salts.

2. BRIEF DESCRIPTION



Main assemblies of the unit are:

- a) furnace assembly;
- b) stand assembly.

The furnace assembly comprises cylindrical shell 4 and cover 3. Inside the shell, on a foamed-chamotte cushion is secured ceramic tube 2 into which is inserted spiral heating element made of 34-626 alloy.


Both the cushion and the ceramic tube are inserted in an asbestos cylinder. The space between the asbestos cylinder and the ceramic tube is filled up with foamed chamotte fines, whereas the space between the asbestos cylinder and the shell is filled up with lagging.

- 3 -

The furnace stand assembly comprises a metal housing 5 the upper cover of which mounts a step-down transformer. The type ААТФ-1 autotransformer is mounted on the housing side wall.

The autotransformer connected into the circuit of the furnace enables to regulate the consumption of energy by the furnace, thus ensuring higher accuracy of temperature control in the furnace working room.

3. SCHEMATIC DIAGRAM OF FURNACE CONNECTION TO SUPPLY



A - Type ААТФ-1 Autotransformer

T - Step-down transformer, 220/22 Volts

H - Heater

4. SPECIFICATION

Supply voltage	220 Volts
Consumption	1850 Watts
Maximum voltage at heater terminals	22 Volts
Maximum working temperature	1200°C
Time required to raise temperature up to 1200°C	120 minutes

- 4 -

Temperature uniformity in the furnace working room	$\pm 5^{\circ}\text{C}$
Stretch of uniform temperature zone	80 mm.
Dimensions of furnace working room:	
Diameter	72 mm.
Height	175 mm.
Overall dimensions of furnace:	
Height	765 mm.
Diameter	455 mm.
Furnace net weight	93 kg.

5. PREPARATION OF FURNACE FOR OPERATION

- a) Unpack the furnace;
- b) Examine the furnace outside and inside the working room;
- c) In case no damage has been found, ground the furnace shell by connecting the terminal "Ground" to the grounding system of the building;
- d) Connect the furnace terminals to the supply mains according to the schematic diagram given in clause 3;
- e) The furnace must be dried out at a temperature ranging within 700 and 900° C for 4 or 5 hours before putting it into service.

6. OPERATION

- a) Set the type AATP-1 autotransformer handle at position "0";
- b) Switch the furnace in;
- c) By gradually turning the autotransformer handle in the direction of voltage increase, check the temperature increase inside the furnace working room using the platinum-platinum-rhodium thermocouple;

- 5 -

d) After the required temperature in the working room has been attained the temperature control must be effected by means of the autotransformer.

e) With the purpose of saving the type J ATP-1 autotransformer from damage it is necessary that, the current consumed from the supply at continuous service of the furnace, would not exceed 8.5 Amps. At this, the current consumed by the furnace will be in the order of 7.5 Amps. Cut in an ammeter with scale up to 10 Amperes to check the current consumed by the furnace;

f) Temperature rise is defined by the speed of heating of the furnace working room, i.e., by the power incoming to the heater.

Best control conditions are ensured by proper selection of the power supplied to the heater, which may be attained with the help of the type J ATP-1 autotransformer.

It is defined experimentally that the temperature of 1200°C is maintained stable when the autotransformer handle is set at the position marked 175 ± 5 Volts.

With ageing of the heater in the course of a long-lasting continuous operation of the furnace, the position of the autotransformer handle for the respective rated temperature may become altered.

7. STORAGE OF FURNACE

The furnace should be stored in dry premises. The ambient air should bear no harmful admixtures resulting in corrosion.

8. DELIVERY VOLUME

a) One Model TP-3 Crucible Type Electric Furnace Set including Type J ATP-1 autotransformer and Step-Down Transformer mounted on the Furnace Stand 1 unit.

- 6 -

b) Technical documents comprising:

Certificate of Works Control Department	1 copy
Operation and maintenance instructions	1 copy

9. SPARE PARTS

a) Heater assembly including alundum bushing	1 set
b) Ceramic cushion support	1 piece
c) Ceramic guide	1 piece
d) Ceramic socket	1 piece
e) Crucible support (alundum)	1 piece

10. GUARANTEE

The Works guarantees normal operation of the furnace for 500 hours after being put into service, provided the rules of the present Instructions are strictly observed and the furnace has been kept in storage, prior to its use, for no more than 4 months under conditions as stipulated in clause 7 of these Instructions. If the storage time has exceeded 4 months, the guarantee term for a trouble-free furnace operation is cut down accordingly.

"PLATINOPRIBOR" STATE WORKS.

FOR OFFICIAL USE ONLY

LIST OF SPARE PARTS
FOR ARTICLE TYPE VI-3
/ INSTEAD OF PART 9/

No	Spare parts	Index of spare parts	Unit	Quantity of spare parts necessary for two years work
1	2	3	4	5
1.	Block	1620106	pos	2
2.	Block	1620107	-"-	2
3.	Substitution	1620108	-"-	2
4.	Ceramic box	1620109	-"-	2
5.	Guide	1620110	-"-	2
6.	Push	1620111	-"-	2
7.	Plug	1620112	-"-	1
8.	Substitution	1620113	-"-	2
9.	Lower packing	1620114	-"-	3
10.	Upper packing	1620115	-"-	2
11.	Ceramic box	1620116	-"-	2
12.	Header 34-626			
	10 5.4	1620117	kg.	1.476

FOR OFFICIAL USE ONLY

25X1

Page Denied